


PL-01026 / PL-01077

## POLYETHYLENE 3-DIMENSIONAL BAGS


ASI's 3-D single-use bags have been specifically designed for larger scale applications. 3-D bags are beneficial to larger scale applications because their designs and shapes conform to vessels better and present fewer creases than traditional 2-D bags. They are made from high quality film that carries a complete testing package, as well as Non-Animal Origin status. These systems are provided completely packaged and sterile, ready for immediate use.






- Port 1: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp
- Port 2: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp
- Port 3: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp

| ASI Catalog # | Size      | Tank #  |
|---------------|-----------|---------|
| SS00046-I     | 50 Liter  | TK-0050 |
| SS00047-I     | 100 Liter | TK-0100 |
| SS00048-I     | 200 Liter | TK-0200 |
| SS00049-I     | 500 Liter | TK-0500 |




- Port 1: 1/2" ID x 3/4" OD Clear C-Flex 374 (60")  
MPX Male + Cap, Pinch Clamp
- Port 2: 1/8" ID x 1/4" OD Clear C-Flex 374 (40")  
Female Luer + Swabable Valve + Cap, Pinch Clamp
- Port 3: 1/2" ID x 3/4" OD Clear C-Flex 374 (60")  
1-1/2" Ladish Tri-Clover, Pinch Clamp
- Port 4: 1/2" ID x 3/4" OD Clear C-Flex 374 (60")  
1-1/2" Ladish Tri-Clover, Pinch Clamp

| ASI Catalog # | Size      | Tank #  |
|---------------|-----------|---------|
| SS00050-I     | 50 Liter  | TK-0050 |
| SS00051-I     | 100 Liter | TK-0100 |
| SS00052-I     | 200 Liter | TK-0200 |
| SS00053-I     | 500 Liter | TK-0500 |



- Port 1: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp
- Port 2: 1/2" ID x 3/4" OD Silicone Tubing (36")  
ReadyMate DAC w/ 1/2" Barb, Pinch Clamp
- Port 3: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp


| ASI Catalog # | Size      | Tank #  |
|---------------|-----------|---------|
| SS00054-I     | 50 Liter  | TK-0050 |
| SS00055-I     | 100 Liter | TK-0100 |
| SS00056-I     | 200 Liter | TK-0200 |
| SS00057-I     | 500 Liter | TK-0500 |



- Port 1: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp
- Port 2: 1/2" ID x 3/4" OD Silicone Tubing (36")  
Female Kleenpak Connector with 1/2" Barb, Pinch Clamp
- Port 3: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp


| ASI Catalog # | Size      | Tank #  |
|---------------|-----------|---------|
| SS00058-I     | 50 Liter  | TK-0050 |
| SS00059-I     | 100 Liter | TK-0100 |
| SS00060-I     | 200 Liter | TK-0200 |
| SS00061-I     | 500 Liter | TK-0500 |

PL-01026 / PL-01077 POLYETHYLENE 3-D BAGS




- Port 1: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp
- Port 2: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp
- Port 3: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp

| ASI Catalog # | Size      | Tank #     |
|---------------|-----------|------------|
| SS00124-I     | 50 Liter  | TK-0050-06 |
| SS00125-I     | 100 Liter | TK-0100-06 |
| SS00126-I     | 200 Liter | TK-0200-06 |
| SS00127-I     | 500 Liter | TK-0500-06 |




- Port 1: 1/2" ID x 3/4" OD Clear C-Flex 374 (60")  
MPX Male + Cap, Pinch Clamp
- Port 2: 1/8" ID x 1/4" OD Clear C-Flex 374 (40")  
Female Luer + Swabable Valve + Cap, Pinch Clamp
- Port 3: 1/2" ID x 3/4" OD Clear C-Flex 374 (60")  
1-1/2" Ladish Tri-Clover, Pinch Clamp
- Port 4: 1/2" ID x 3/4" OD Clear C-Flex 374 (60")  
1-1/2" Ladish Tri-Clover, Pinch Clamp

| ASI Catalog # | Size      | Tank #     |
|---------------|-----------|------------|
| SS00128-I     | 50 Liter  | TK-0050-06 |
| SS00129-I     | 100 Liter | TK-0100-06 |
| SS00130-I     | 200 Liter | TK-0200-06 |
| SS00131-I     | 500 Liter | TK-0500-06 |




- Port 1: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp
- Port 2: 1/2" ID x 3/4" OD Silicone Tubing (36")  
ReadyMate DAC with 1/2" Barb, Pinch Clamp
- Port 3: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp

| ASI Catalog # | Size      | Tank #     |
|---------------|-----------|------------|
| SS00132-I     | 50 Liter  | TK-0050-06 |
| SS00133-I     | 100 Liter | TK-0100-06 |
| SS00134-I     | 200 Liter | TK-0200-06 |
| SS00135-I     | 500 Liter | TK-0500-06 |




- Port 1: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp
- Port 2: 1/2" ID x 3/4" OD Silicone Tubing (36")  
Female Kleenpak Connector with 1/2" Barb, Pinch Clamp
- Port 3: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp

| ASI Catalog # | Size      | Tank #     |
|---------------|-----------|------------|
| SS00136-I     | 50 Liter  | TK-0050-06 |
| SS00137-I     | 100 Liter | TK-0100-06 |
| SS00138-I     | 200 Liter | TK-0200-06 |
| SS00139-I     | 500 Liter | TK-0500-06 |




- Port 1: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp
- Port 2: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp
- Port 3: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp

| ASI Catalog # | Size       | Tank #     |
|---------------|------------|------------|
| SS00140-I     | 200 Liter  | TK-0250-08 |
| SS00141-I     | 500 Liter  | TK-0500-08 |
| SS00142-I     | 1000 Liter | TK-1000-08 |




- Port 1: 1/2" ID x 3/4" OD Clear C-Flex 374 (60")  
MPX Male + Cap, Pinch Clamp
- Port 2: 1/8" ID x 1/4" OD Clear C-Flex 374 (40")  
Female Luer + Swabable Valve + Cap, Pinch Clamp
- Port 3: 1/2" ID x 3/4" OD Clear C-Flex 374 (60")  
1-1/2" Ladish Tri-Clover, Pinch Clamp
- Port 4: 1/2" ID x 3/4" OD Clear C-Flex 374 (60")  
1-1/2" Ladish Tri-Clover, Pinch Clamp

| ASI Catalog # | Size       | Tank #     |
|---------------|------------|------------|
| SS00143-I     | 200 Liter  | TK-0250-08 |
| SS00144-I     | 500 Liter  | TK-0500-08 |
| SS00145-I     | 1000 Liter | TK-1000-08 |



- Port 1: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp
- Port 2: 1/2" ID x 3/4" OD Silicone Tubing (36")  
ReadyMate DAC with 1/2" Barb, Pinch Clamp
- Port 3: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp

| ASI Catalog # | Size       | Tank #     |
|---------------|------------|------------|
| SS00146-I     | 200 Liter  | TK-0250-08 |
| SS00147-I     | 500 Liter  | TK-0500-08 |
| SS00148-I     | 1000 Liter | TK-1000-08 |



- Port 1: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp
- Port 2: 1/2" ID x 3/4" OD Silicone Tubing (36")  
Female Kleenpak Connector with 1/2" Barb, Pinch Clamp
- Port 3: 1/2" ID x 3/4" OD Silicone Tubing (20")  
MPX Male + Cap, Pinch Clamp

| ASI Catalog # | Size       | Tank #     |
|---------------|------------|------------|
| SS00149-I     | 200 Liter  | TK-0250-08 |
| SS00150-I     | 500 Liter  | TK-0500-08 |
| SS00151-I     | 1000 Liter | TK-1000-08 |